# Installation Manual

## CuratOR<sup>®</sup> EX4342 Color LCD Monitor

### Important

Please read this "Installation Manual", and "Instructions for Use" (separate volume) carefully to familiarize yourself with safe and effective usage.

- 
- Refer to the "Instructions for Use" for information on the installation / connection of the monitor. The "Installation Manual" is part of the "Instructions for Use".
-

This product has been adjusted specifically for use in the region to which it was originally shipped. If operated outside this region, the product may not perform as stated in the specifications.

---

No part of this manual may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, or otherwise, without the prior written permission of EIZO GmbH. EIZO GmbH is under no obligation to hold any submitted material or information confidential unless prior arrangements are made pursuant to EIZO GmbH's receipt of said information. Although every effort has been made to ensure that this manual provides up-to-date information, please note that EIZO product specifications are subject to change without notice.

---

# CONTENTS

<b>CONTENTS</b> .....	<b>3</b>
<b>Chapter 1 Basic Adjustment / Setting</b> .....	<b>4</b>
1-1. Loading the Monitor Settings.....	4
● Saving current settings to “User1” through “User5” .....	5
1-2. Switching Input Signals .....	5
1-3. Using Custom Keys .....	6
● Setting a function to custom keys .....	6
1-4. Locking Operation Buttons .....	6
<b>Chapter 2 Advanced Adjustments / Settings</b> ...7	
2-1. Basic Operation of the Setting Menu .....	7
2-2. Setting Menu Functions .....	9
● Color.....	9
● Signal .....	10
● Preferences.....	11
● Information .....	13
● PbyP.....	13
● PinP.....	14
<b>Chapter 3 Troubleshooting</b> .....	<b>15</b>
3-1. No Picture Displayed .....	15
3-2. Imaging Problems.....	16
3-3. Other Problems .....	16
<b>Chapter 4 Reference</b> .....	<b>17</b>
4-1. Main Default Settings .....	17
4-2. Pin Specifications.....	18
<b>Appendix</b> .....	<b>19</b>
Trademark .....	19
License .....	19
LIMITED WARRANTY.....	20
Recycling Information .....	21

# Chapter 1 Basic Adjustment / Setting

This chapter describes the basic functions that can be adjusted and set by pressing the buttons on the front of the monitor.

## 1-1. Loading the Monitor Settings


Basic settings are made in advance at the factory.

It is possible to switch to the following preset settings that have been set in advance.

DisplayPort [sRGB] DisplayPort [BT.2020]	Basic settings for DisplayPort. The input signal is set to DisplayPort when this item is selected. Settings cannot be changed. For information about the main default settings, see <a href="#">“4-1. Main Default Settings” (page 3)</a> .
HDMI [BT.2020]	Basic settings for HDMI. The input signal is set to HDMI when this item is selected. Settings cannot be changed. For information about the main default settings, see <a href="#">“4-1. Main Default Settings” (page 3)</a> .
SDI [BT.2020]	Basic settings for SDI 1. The input signal is set to SDI 1 when this item is selected. Settings cannot be changed. For information about the main default settings, see <a href="#">“4-1. Main Default Settings” (page 3)</a> .
“User1” through “User5”	You can save the current setting values and display status in the following items. These items are not selectable unless they have been saved. For details, see <a href="#">“Saving current settings to “User1” through “User5” (page 4)</a> .

### 1. Press PRESET.

The preset menu is displayed.

### 2. Press F1/Left, F2/Right, F3/Down and F4/Up to select the item to load, and then press ENTER.

The pre-made settings are loaded and the settings are changed.

#### Note

- The preset menu will disappear after five seconds have elapsed if no operation is performed.

## ● Saving current settings to “User1” through ”User5”

You can save the current setting values and display status to “User1” through ”User5”.

### 1. Press MENU.

The Setting menu appears.

### 2. Press F3/Down or F4/Up to select “Preferences”, and then press ENTER.

The “Preferences” menu appears.

### 3. Press F3/Down or F4/Up to select “PRESET Save”, and then press ENTER.

The “PRESET Save” menu appears.

### 4. Press F3/Down or F4/Up to select one from “User1” to “User5”, and then press ENTER.


Current setting values and display status will be saved in the selected item.

## 1-2. Switching Input Signals

Select the signal input to the monitor and display it on the screen.

### 1. Press INPUT.

The input port menu appears.

### 2. Press F3/Down or F4/Up to select an input signal.

This switches to the selected input signal.

### 3. Press INPUT.

This exits the menu screen.

---

#### Note

- If you press INPUT while the input port menu is displayed in dual screen mode (PbyP or PinP), the “PbyP” / “PinP” menu is displayed. For details, see [“PbyP” \(page 9\)](#) and [“PinP” \(page 10\)](#).
  - After the input port menu is displayed, the input port menu will disappear after 10 seconds if no key operations are performed.
-

## 1-3. Using Custom Keys

---

Execute the functions assigned to the custom keys (F1/Left, F2/Right, F3/Down and F4/Up). For information on how to set a function to each custom key, see [“Setting a function to custom keys”](#) (page 5).

### 1. Press F1/Left, F2/Right, F3/Down and F4/Up.

The function associated with the custom key that is pressed will run.

---

**Note**

- For dual screen display (PbyP or PinP), adjust the “Active Window” side.
- 

### ● Setting a function to custom keys

#### 1. Press MENU.

The Setting menu appears.

#### 2. Press F3/Down or F4/Up to select “Preferences”, and then press ENTER.

The “Preferences” menu appears.

#### 3. Press F3/Down or F4/Up to select “Custom Key”, and then press ENTER.

The “Custom Key” menu appears.

#### 4. Press F3/Down or F4/Up to select which custom key to set the function to, and then press ENTER.

#### 5. Press F3/Down or F4/Up to select the function, and then press ENTER.

The selected function is assigned to the corresponding custom key.

## 1-4. Locking Operation Buttons

---

Adjusted or set states can be made unchangeable.

### Procedure

1. Press the MENU and ENTER keys for at least 4 seconds.

---

**Note**

- Press the MENU and ENTER keys for at least 4 seconds to exit Key Lock.
-


# Chapter 2 Advanced Adjustments / Settings

This chapter describes advanced monitor adjustments and settings procedures using the main menu.


## 2-1. Basic Operation of the Setting Menu

### 1. Menu display

1. Press MENU.

The Setting menu appears.


## 2. Adjustments / Settings

Pressing MENU during adjustment/setting will cancel the adjustment/setting and restore the state prior to making changes.


1. Press F3/Down or F4/Up to select a menu to adjust/set, and then press ENTER.

The submenu appears.

2. Press F3/Down or F4/Up to select an item to adjust/set, and then press ENTER.

The Adjustment/Setting menu appears.

3. Press F1/Left or F2/Right to adjust/set the selected item, and then press ENTER to confirm the settings.

The submenu appears.

## 3. Exiting

1. Press MENU.

The Setting menu appears.

2. Press MENU.

Exits the Setting menu.

## 2-2. Setting Menu Functions


### ● Color

Color (All)		
Active Window	[	All ]
Backlight	[	100 ]
Brightness	[	0 ]
Contrast	[	100 ]
Temperature	[	6500K ]
Gamma	[	2.4 ]
Color Gamut	[	sRGB ]
Hue	[	0 ]
Saturation	[	0 ]
Gain	[	0 ]

Function	Settings	Description
Active Window* <sup>1</sup>	PbyP Left Window / Right Window	Select which window to adjust settings when dual screen display (PbyP or PinP).
	PinP Large Window / Small Window	<b>Note</b> • For dual screen display (PbyP or PinP), color adjustment can be performed for each screen selected in “Active Window”.
	Single All	
Backlight	0 to 100	Adjust screen brightness by changing the brightness of the backlight (light source from the LCD back panel). <b>Note</b> • For dual screen display (PbyP or PinP), settings are applied to “Left Window” / “Large Window”. • During PbyP/PinP, the backlight is disabled when either “Right Window” or “Small Window” are selected in the Active Window.
Brightness	0 to 100	Adjust the brightness of black.
Contrast	0 to 100	Contrast can be adjusted. Adjust the brightness of the screen by changing the video signal level.
Temperature	Native / 6500K / 9300K / User	Adjust the color temperature. The color temperature is normally used to express the hue of “White” and / or “Black” with a numerical value. The value is expressed in degrees “K” (Kelvin). The screen becomes reddish at low color temperatures and bluish at high color temperatures, similar to the temperatures of a flame. <b>Note</b> • Select “Native” to display the original colors of the LCD panel. • “Gain” allows you to perform more advanced adjustment. When gain is changed, “Temperature” is changed to “User”. • The gain preset values are set for each color temperature setting value.
Gamma	1.8 / 2.0 / 2.2 / 2.4 / 2.6 / DICOM / HLG / PQ	Adjust the gamma. While the brightness of the monitor varies depending on the input signal, the variation rate is not proportional to the input signal. Maintaining the balance between the input signal and brightness of the monitor is referred to as “Gamma correction”.
Color Gamut	Native / sRGB / BT.2020	Set the color reproduction area (color gamut). Color Gamut is the range of colors that can be represented. Multiple standards are defined. <b>Note</b> • Select “Native” to display the original color gamut of the LCD panel.
Hue	-100 to 100	Adjust the hue.
Saturation	-100 to 100	Adjust the color saturation.
Gain	0 to 255	The brightness of each color component (red, green, and blue) is called “gain”. The hue of “white” can be changed by adjusting the gain.

\*1 Disabled when single screen display.

## ● Signal




Function	Settings	Description
Active Window* <sup>1</sup>	(When single screen display mode.) All (When "Multi Window" is set to "PbyP".) Left Window / Right Window (When "Multi Window" is set to "PinP".) Large Window / Small Window	Select which window to adjust settings when dual screen display (PbyP or PinP).
Input Color Format* <sup>2</sup>	(When input signal is DisplayPort.) Auto / YUV / RGB (When input signal is HDMI.) Auto / YUV 4:2:2 / YUV 4:4:4 / RGB (When input signal is SDI.) Auto / YCC4:2:2 10bit / YCC4:2:2 12bit / YCC4:4:4 10bit / YCC4:4:4 12bit / RGB4:4:4 10bit / RGB4:4:4 12bit	The color space of the input signal can be specified. Try changing this setting if colors are not displayed correctly. Normally, select "Auto". When an SDI signal is input, the bit depth is specified using the "Input Color Format". If the specified bit depth is invalid, a signal error is displayed.
YUV Color Matrix* <sup>2</sup>	Auto / BT.601 / BT.709 / BT.2020	Select the YUV format of the input signal. Use this setting for gradation collapse or other occurrences where the screen is displayed incorrectly because of a problem with the input signal.
Input Range	(When input signal is DisplayPort or HDMI.) Auto / Full / Limited (When input signal is SDI or DVI.) Full / Limited	Depending on the video reproduction device, there may be a restriction on black and white video signal levels output to the monitor. This kind of signal is called "Limited range". On the other hand, unlimited signals are called "Full range". <ul style="list-style-type: none"> <li>• "Auto" Input signal brightness range is automatically judged and is displayed appropriately (recommended setting). Depending on the video reproduction device, Limited range and Full range may not be able to be judged by the monitor. In such a case, selecting "Full" or "Limited" will allow it to be displayed appropriately.</li> <li>• "Full" To be selected in case of Full range signals. Appropriate display can be obtained when this is selected, in case both blacks and whites are corrupted.</li> <li>• "Limited" To be selected in case of Limited range signals. When this is selected, the output signal range is expanded from 0 to 255 to obtain an appropriate display when black is pale and white is dull.</li> </ul>

\*1 Disabled when single screen display.

\*2 Disabled when input signal is DVI.

## ● Preferences

Function	Settings	Description						
Multi Window* <sup>1</sup>	Single / PbyP / PinP	<p>Select single screen display or dual screen display (PbyP or PinP).</p> <ul style="list-style-type: none"> <li>• “Single” Display selected input signal on the screen.</li> <li>• “PbyP” Display two different signals side-by-side on one screen simultaneously.</li> <li>• “PinP” Display two different signals on one screen by overlapping a subscreen on a main screen.</li> </ul> <p><b>Note</b></p> <ul style="list-style-type: none"> <li>• In PbyP and PinP, the same input signal can be displayed on two screens. Example: During PbyP display, select “DisplayPort” for both “Left Window” and “Right Window” on the “PbyP” screen.</li> <li>• “SDI 1” and “SDI 2” cannot be displayed simultaneously.</li> </ul>						
Display Direction* <sup>2</sup>	Normal / 180° / Mirror	<p>Select the rotation angle or horizontally flip the main screen.</p> <p><b>Note</b></p> <ul style="list-style-type: none"> <li>• Displays  on the screen when “180°” selected.</li> <li>• Displays  on the screen when “Mirror” selected.</li> </ul>						
Languages	English / German / French / Spanish / Italian / Swedish / Japanese / Simplified Chinese / Traditional Chinese / Russian / Portuguese	Select the language to be used for display in settings and all menus.						
PRESET Save	<table border="1"> <tr> <td>User1</td> <td rowspan="5">OK / Cancel</td> </tr> <tr> <td>User2</td> </tr> <tr> <td>User3</td> </tr> <tr> <td>User4</td> </tr> <tr> <td>User5</td> </tr> </table>	User1	OK / Cancel	User2	User3	User4	User5	<p>Save current setting values and display status for each user.</p> <p><b>Note</b></p> <ul style="list-style-type: none"> <li>• Setting items identical to preset settings (DisplayPort [sRGB]/DisplayPort [BT.2020]/HDMI [BT.2020]/SDI [BT.2020]) are saved. For the setting items of preset settings, see “4-1. Main Default Settings” (page 17).</li> </ul>
User1	OK / Cancel							
User2								
User3								
User4								
User5								
Custom Key	<table border="1"> <tr> <td>[F1]</td> <td rowspan="4">Off / Backlight / Gamma / Color Gamut / Multi Window / Display Direction / Swap</td> </tr> <tr> <td>[F2]</td> </tr> <tr> <td>[F3]</td> </tr> <tr> <td>[F4]</td> </tr> </table>	[F1]	Off / Backlight / Gamma / Color Gamut / Multi Window / Display Direction / Swap	[F2]	[F3]	[F4]	Select the function that is assigned to each custom key. For details, see “Setting a function to custom keys” (page 6).	
[F1]	Off / Backlight / Gamma / Color Gamut / Multi Window / Display Direction / Swap							
[F2]								
[F3]								
[F4]								

Function		Settings		Description	
Advanced Settings	Auto Input Detection <sup>*2</sup>	On / Off		<p>This function automatically recognizes the connector through which computer signals are input, and displays images on the screen accordingly.</p> <ul style="list-style-type: none"> <li>• “On” When the monitor is connected to multiple computers, if a specific computer enters power saving mode or no signals are input to the monitor, the connector is automatically changed to another one to which signals are input.</li> <li>• “Off” Set to this when manually selecting input signals. In this case, you can select the input signal to display using the INPUT. For details, see <a href="#">“1-2. Switching Input Signals” (page 5)</a>.</li> </ul> <p><b>Note</b></p> <ul style="list-style-type: none"> <li>• When this function is set to “On”, the monitor only enters the power saving mode when signals are not input from any computer.</li> </ul>	
	Input Skip <sup>*2</sup>	DisplayPort	- / Skip		<p>When “Auto Input Detection” is “On”, select the port to skip signal detection.</p> <p><b>Note</b></p> <ul style="list-style-type: none"> <li>• If all of the ports are set to “Skip”, the operation will be the same as when “Auto Input Detection” is “Off”. Detection of input signals is not performed.</li> </ul>
		HDMI			
		SDI 1			
		SDI 2			
		DVI			
On-Screen Logo	On / Off		Select whether to display or hide the logo when the monitor starts.		
SDI1 Input Equalizer	On / Off		<p>Emphasizes the input signal of SDI 1.</p> <p><b>Note</b></p> <ul style="list-style-type: none"> <li>• Set to “Off” if noise occurs using the SDI short cable.</li> </ul>		
RS-232C	19200 bps / 115200 bps		Specifies the transmission rate.		
Monitor Reset	OK / Cancel		Restores all settings to their default values.		

\*1 Disabled when “Display Direction” is set to “180°” or “Mirror”.

\*2 Disabled when dual screen display (PbyP or PinP).

## ● Information

Display the monitor information (model name, serial number (S/N), usage time) and the input signal information.

**Example:**

```

Information
CuratOR EX4342      S/N: 12345678
Usage Time (h)      0

DisplayPort
1920 X 1080
fH: 67.50 kHz
fV: 60.00 Hz
fD: 148.5 MHz
  
```

## ● PbyP

If you press INPUT while in PbyP dual screen mode, the “PbyP” menu is displayed.


```

PbyP Settings
Left Window [ DVI ]
Right Window [ SDI 1 ]
Swap
  
```

Function	Setting	Description
Left Window	DisplayPort / HDMI / SDI 1 / SDI 2 / DVI	Select the signal input to the left window when PbyP dual screen display. <b>Note</b> • “SDI 1” and “SDI 2” cannot be displayed simultaneously.
Right Window	DisplayPort / HDMI / SDI 1 / SDI 2 / DVI	Select the signal input to the right window when PbyP dual screen display. <b>Note</b> • “SDI 1” and “SDI 2” cannot be displayed simultaneously.
Swap		Swap the input ports of left window and right window.

## ● PinP


If you press INPUT while in PinP dual screen mode, the “PinP” menu is displayed.

Function	Setting	Description
Large Window	DisplayPort / HDMI / SDI 1 / SDI 2 / DVI	Select the signal input to the large window when PinP dual screen display. <b>Note</b> • “SDI 1” and “SDI 2” cannot be displayed simultaneously.
Small Window	DisplayPort / HDMI / SDI 1 / SDI 2 / DVI	Select the signal input to the small window when PinP dual screen display. <b>Note</b> • “SDI 1” and “SDI 2” cannot be displayed simultaneously.
PinP Size	Small / Large	Select size of the small window.
PinP Position	Upper Right / Lower Right / Lower Left / Upper Left	Select position of the small window.
Swap		Swap the input ports of large window and small window.

# Chapter 3 Troubleshooting

## 3-1. No Picture Displayed

Problem	Possible cause and remedy
<p><b>1. No picture.</b></p>	<ul style="list-style-type: none"><li>• Check whether the power cord is connected properly.</li><li>• Check whether the DC OUT terminal and the DC IN terminal are connected properly.</li><li>• Turn on the power switch.</li><li>• Check whether the main power indicator of the AC adapter is on.</li><li>• Turn off the power and then turn it on again.</li></ul>
<p><b>2. The message below appears.</b></p> <ul style="list-style-type: none"><li>• This message appears when no signal is input. Example: </li><li>• This message indicates that the input signal is outside the specified frequency range. Example: </li></ul>	<p>These messages appear when the signal is not input correctly even though the monitor is functioning properly.</p> <ul style="list-style-type: none"><li>• The message shown on the left may appear, because some devices to be connected do not output the signal immediately after power-on.</li><li>• Check whether the device to be connected is turned on.</li><li>• Check whether the signal cable is connected properly.</li><li>• Turn off the power and then turn it on again.</li></ul> <ul style="list-style-type: none"><li>• Check whether the device to be connected is configured to meet the resolution and vertical scan frequency requirements of the monitor.</li><li>• Reboot the device to be connected.</li></ul>

## 3-2. Imaging Problems

Problem	Possible cause and remedy
1. The screen is too bright or too dark.	<ul style="list-style-type: none"> <li>Use "Brightness" or "Backlight" in the Setting menu to adjust it.</li> </ul>
2. Text is blurred	<ul style="list-style-type: none"> <li>Check whether the computer is configured to meet the resolution and vertical scan frequency requirements of the monitor (see "Displayable Input Signals" in the Instructions for Use).</li> </ul>
3. Afterimages appear	<ul style="list-style-type: none"> <li>Afterimages are particular to LCD monitors. Avoid displaying the same image for a long time.</li> <li>Use the screen saver or power saving function to avoid displaying the same image for extended periods of time.</li> <li>An afterimage may appear even after a short time period has elapsed depending on the displayed image. If this occurs, changing the image or leaving the power off for a few hours may solve the problem.</li> </ul>
4. Green, red, blue, white dots or dark dots remain on the screen	<ul style="list-style-type: none"> <li>This is due to LCD panel characteristics and not a malfunction.</li> </ul>
5. The screen suddenly went black. / The screen does not become brighter after increasing the "Backlight" value in the Settings menu.	<ul style="list-style-type: none"> <li>This product is equipped with a heat dissipation fan. If the fan is not operating properly, screen brightness can be limited in order to lower the internal temperature. In this case, please make sure that the fan is operating properly.</li> </ul>

## 3-3. Other Problems

Problem	Possible cause and remedy
1. The Setting menu does not appear.	<ul style="list-style-type: none"> <li>Make sure the lock function is not set to "On." (See "1-4. Locking Operation Buttons" (page 5).)</li> </ul>
2. The message below appears. <ul style="list-style-type: none"> <li>This will be displayed when a 6-bit (64 gradient) signal is being input. Example:</li> </ul> <div data-bbox="204 1301 708 1503" style="border: 1px solid black; background-color: #333; color: #eee; padding: 5px; margin-top: 10px;"> <p style="text-align: center;">Grayscale Display Error</p> <p style="text-align: center;">6-bit signal has been detected. Turn off the main power on the back of the monitor, and then turn it on again.</p> </div>	This will be displayed when the signal is not displaying properly. <ul style="list-style-type: none"> <li>Turn off the power to the connected device and this product and then turn on the power again.</li> <li>Check the settings on the connected device.</li> </ul>

# Chapter 4 Reference

## 4-1. Main Default Settings

### Preset


Item		Default			
		DisplayPort [sRGB]	DisplayPort [BT.2020]	HDMI [BT.2020]	SDI [BT.2020]
Input Port	MAIN (Left Window/Large Window)	DisplayPort	DisplayPort	HDMI	SDI 1
	SUB (Right Window/Small Window)	SDI 1	SDI 1	SDI 1	DVI
Preferences	Multi Window	Single	Single	Single	Single
	Display Direction	Normal	Normal	Normal	Normal
Color	Backlight	100	100	100	100
	Brightness	0	0	0	0
	Contrast	100	100	100	100
	Temperature	6500K	6500K	6500K	6500K
	Gamma	2.2	2.2	2.2	2.2
	Color Gamut	sRGB	BT.2020	BT.2020	BT.2020
	Saturation	0	0	0	0
Signal	Hue	0	0	0	0
	Input Color Format	Auto	Auto	Auto	Auto
	YUV Color Matrix	Auto	BT.2020	BT.2020	BT.2020
	Input Range	Auto	Auto	Auto	Limited

### Others

Preferences	Language	English
	On-Screen Logo	On
Custom Key	F1	Off
	F2	Off
	F3	Off
	F4	Off
Key Lock		Unlock


## 4-2. Pin Specifications

### ● DC IN terminal

Pin No.	Signal name
1	GND
2	GND
3	F.E
4	48 VDC
5	48 VDC

### ● DC OUT terminal on the main unit


Pin No.	Signal name
1	GND
2	GND
3	5 VDC
4	5 VDC

HIROSE P/N: HR10A-7R-4S(74)

Recommended plug

HIROSE P/N: HR10A-7P-4P(74)

### ● DC OUT terminal of AC adapter

Pin No.	Signal name
1	GND
2	GND
3	F.E
4	48 VDC
5	48 VDC

# Appendix

## Trademark

The terms HDMI and HDMI High-Definition Multimedia Interface, and the HDMI Logo are trademarks or registered trademarks of HDMI Licensing, LLC in the United States and other countries.

The DisplayPort Compliance Logo and VESA are registered trademarks of the Video Electronics Standards Association.

The SuperSpeed USB Trident Logo is a registered trademark of USB Implementers Forum, Inc.

The USB Power Delivery Trident Logos are trademarks of USB Implementers Forum, Inc.

USB Type-C and USB-C are registered trademarks of USB Implementers Forum, Inc.

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communications of medical information.

Kensington and Microsaver are registered trademarks of ACCO Brands Corporation.

Thunderbolt is a trademark of Intel Corporation in the United States and/or other countries.

Microsoft and Windows are registered trademarks of Microsoft Corporation in the United States and other countries.

Adobe is a registered trademark of Adobe Systems Incorporated in the United States and other countries.

Apple, macOS, Mac OS, OS X, Macintosh and ColorSync are registered trademarks of Apple Inc.

ENERGY STAR is a registered trademark of the United States Environmental Protection Agency in the United States and other countries.

EIZO, the EIZO Logo, ColorEdge, CuratOR, DuraVision, FlexScan, FORIS, RadiCS, RadiForce, RadiNET, Raptor and ScreenManager are registered trademarks of EIZO Corporation in Japan and other countries.

ColorEdge Tablet Controller, ColorNavigator, EcoView NET, EIZO EasyPIX, EIZO Monitor Configurator, EIZO ScreenSlicer, G-Ignition, i-Sound, Quick Color Match, RadiLight, Re/Vue, SafeGuard, Screen Administrator, Screen InStyle, ScreenCleaner and UniColor Pro are trademarks of EIZO Corporation.

All other company names, product names, and logos are trademarks or registered trademarks of their respective owners.

## License

The bitmap font used for this product is designed by Ricoh Industrial Solutions Inc.

# LIMITED WARRANTY

EIZO GmbH (hereinafter referred to as “**EIZO**”) and distributors authorized by EIZO (hereinafter referred to as the “**Distributors**”) warrant, subject to and in accordance with the terms of this limited warranty (hereinafter referred to as the “**Warranty**”), to the original purchaser (hereinafter referred to as the “**Original Purchaser**”) who purchased the product specified in this document (hereinafter referred to as the “**Product**”) from EIZO or Distributors, that EIZO and Distributors shall, at their sole discretion, either repair or replace the Product at no charge if the Original Purchaser becomes aware within the Warranty Period (defined below) that the Product malfunctions or is damaged in the course of normal use of the Product in accordance with the description in the instruction manual of the Product (hereinafter referred to as the “**User’s Manual**”).

EIZO and Distributors shall bear no liability or obligation with regard to the Product in relation to the Original Purchaser or any third parties other than as provided under this Warranty.

In repairing the Product, EIZO and Distributors will use renewal parts which comply with our QC standards. If the unit cannot be repaired due to its condition or the stockout of a relevant part, EIZO and Distributors may offer the replacement by a product with equivalent performance instead of repairing it.

The Warranty is valid only in the countries or territories where the Distributors are located. The Warranty does not restrict any legal rights of the Original Purchaser.

Notwithstanding any other provision of this Warranty, EIZO and Distributors shall have no obligation under this Warranty whatsoever in any of the cases as set forth below:

- (a) Any defect of the Product caused by freight damage, modification, alteration, abuse, misuse, accident, incorrect installation, disaster, adherent dust, faulty maintenance and/or improper repair by third party other than EIZO and Distributors;
- (b) Any incompatibility of the Product due to possible technical innovations and/or regulations;
- (c) Any deterioration of the sensor, including measurement value of the sensor;
- (d) Any defect of the Product caused by external equipments;
- (e) Any defect of the Product caused by the use under environmental conditions other than EIZO and Distributors’ estimation;
- (f) Any deterioration of the attachments of the Product (e.g. cables, User’s Manual, CD-ROM, etc.);
- (g) Any deterioration of the consumables, and/or accessories of the Product (e.g. batteries, remote controller, touch pen, etc.);
- (h) Any deformation, discoloration, and/or warp of the exterior of the Product including that of the surface of the LCD panel;
- (i) Any defect of the Product caused by placement in a location where it may be affected by strong vibration or shocks;
- (j) Any defect of the Product caused by liquid leaking from battery;
- (k) Any deterioration of display performance caused by the deterioration of expendable parts such as the LCD panel and/or backlight, etc. (e.g. changes in brightness, changes in brightness uniformity, changes in color, changes in color uniformity, defects in pixels including burnt pixels, etc.);
- (l) Any deterioration or malfunction of the cooling fan caused by adherent dust.

To obtain service under the Warranty, the Original Purchaser must deliver the Product, freight prepaid, in its original package or other adequate package affording an equal degree of protection, assuming the risk of damage and/or loss in transit, to the local Distributor. The Original Purchaser must present proof of purchase of the Product and the date of such purchase when requesting services under the Warranty.

The Warranty Period for any replaced and/or repaired product under this Warranty shall expire at the end of the original Warranty Period.

EIZO OR DISTRIBUTORS ARE NOT RESPONSIBLE FOR ANY DAMAGE TO, OR LOSS OF, DATA OR OTHER INFORMATION STORED IN ANY MEDIA OR ANY PART OF ANY PRODUCT RETURNED TO EIZO OR DISTRIBUTORS FOR REPAIR.

EIZO AND DISTRIBUTORS MAKE NO FURTHER WARRANTIES, EXPRESSED OR IMPLIED, WITH RESPECT TO, INCLUDING, WITHOUT LIMITATION, THE PRODUCT AND ITS QUALITY, PERFORMANCE, MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR USE. IN NO EVENT SHALL EIZO OR DISTRIBUTORS BE LIABLE FOR ANY INCIDENTAL, INDIRECT, SPECIAL, CONSEQUENTIAL OR OTHER DAMAGE WHATSOEVER (INCLUDING, WITHOUT LIMITATION, DAMAGES FOR LOSS OF PROFIT, BUSINESS INTERRUPTION, LOSS OF BUSINESS INFORMATION, OR ANY OTHER PECUNIARY LOSS) ARISING OUT OF THE USE OR INABILITY TO USE THE PRODUCT OR IN ANY CONNECTION WITH THE PRODUCT, WHETHER BASED ON CONTRACT, TORT, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE, EVEN IF EIZO OR DISTRIBUTORS HAVE BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. THIS EXCLUSION ALSO INCLUDES ANY LIABILITY WHICH MAY ARISE OUT OF THIRD PARTY CLAIMS AGAINST THE ORIGINAL PURCHASER. THE ESSENCE OF THIS PROVISION IS TO LIMIT THE POTENTIAL LIABILITY OF EIZO AND DISTRIBUTORS ARISING OUT OF THIS LIMITED WARRANTY AND/OR SALES.

# Recycling Information

## Recycling Information

This product, when disposed of, is supposed to be collected and recycled according to your country's legislation to reduce environmental burden. When you dispose of this product, please contact a distributor or an affiliate in your country.

The contact addressees are listed on the EIZO website below.

[www.eizoglobal.com](http://www.eizoglobal.com)

For recycling information for customers in Switzerland, please refer to the following website.

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